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***Contrel***

*Controlled Release for the Enhancement of Quality of Life of Women*

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### **II. Composition**

Qualitative composition: Cuprum 99,99% (high purity copper)

Quantitative composition: 280 mg (GyneFix<sup>®</sup> 200); 420 mg (GyneFix<sup>®</sup> 330)

### **III. Description**

GyneFix<sup>®</sup> is an innovative concept in intrauterine device (IUD) technology. It is fixed to the uterine fundus. It's frameless and entirely flexible. GyneFix<sup>®</sup> consists of 4 copper tubes (GyneFix<sup>®</sup> 200) or 6 copper tubes (GyneFix<sup>®</sup> 330), threaded on a length of polypropylene suture material. The tubes are prevented from sliding off by the upper and lower tubes, being crimped onto the thread. The proximal end of the thread is provided with a knot which at insertion is placed in the fundal myometrium with an inserter for anchoring the device. Immediately below the anchoring knot a small stainless steel tube is fixed onto the suture thread. This metal element allows clear visualization on ultrasound examination. The total surface of copper, including the inner and outer surfaces, is 200 mm<sup>2</sup> for the GyneFix<sup>®</sup> 200 and 330 mm<sup>2</sup> for the GyneFix<sup>®</sup> 330. GyneFix<sup>®</sup> is a sterile medical device, intended for single use.

### **IV. Clinical Pharmacology**

Available data indicate that the contraceptive effectiveness of copper IUDs is enhanced by a minute quantity of copper being released continuously into the uterine cavity. The exact mechanism by which copper enhances the contraceptive effect of an IUD has not been conclusively demonstrated. Various hypotheses have been advanced; the major one being that copper placed in the uterus interferes with enzymatic or other processes that regulate blastocyst implantation. Recent studies suggest that copper may play a role in reducing sperm transport within the uterine environment.

### **V. Indication and Clinical Use**

#### **Indication Clinical Use**

GyneFix<sup>®</sup> is indicated for contraception.

Every healthy, non pregnant woman, including nulliparous and adolescent women, can use GyneFix<sup>®</sup>. GyneFix<sup>®</sup> is particularly suitable for women with small uterine cavities. It is recommended to use the GyneFix<sup>®</sup> 200 in women with uterine sound length of 8 cm or less and the GyneFix<sup>®</sup> 330 in women with sound length greater than 8 cm.

#### **Time of insertion**

GyneFix<sup>®</sup> can be inserted at any time during the menstrual cycle. It is evident, however, that the physician should rule out the presence of a pregnancy if insertion of the device is planned in the second half of the menstrual cycle.

GyneFix<sup>®</sup> can be inserted immediately after spontaneous or induced first trimester abortion of less than 10 weeks gestational age.

GyneFix<sup>®</sup> is unsuitable for use in the immediate post placental period. A period of at least 6 weeks should be allowed until the uterus has regained its normal size or the subject has had at least one normal menstrual period. In breastfeeding women, however,

a period of 12 weeks-following full term delivery is recommended because of the higher risk of perforation reported in some studies with traditional IUDs.

## **VI. Use-Effectiveness Studies**

GyneFix<sup>®</sup> has been studied since 1985 and over 15,000 woman-years of experience (international, multicenter, non-comparative and comparative clinical trials, including a large proportion of nulligravid/nulliparous women) have been collected.

### **1. Effectiveness of GyneFix<sup>®</sup>**

GyneFix<sup>®</sup> is a very effective contraceptive device. The pregnancy rates are comparable or lower when compared with the most effective currently used high-load copper devices, and to oral contraceptives and sterilization.

#### ***Use in young parous and nulligravid/nulliparous women***

Studies have shown that young parous and nulligravid/nulliparous women using GyneFix<sup>®</sup> are no more susceptible to pregnancy than other age groups, the annual pregnancy rates being close to zero. The low pregnancy rate with GyneFix<sup>®</sup> is attributed to the high delivery of the copper ions in the uterine cavity.

### **2. Bleeding/pain with GyneFix<sup>®</sup>**

Due to its design characteristics, flexibility, absence of a frame and the reduced space it occupies within the uterine cavity, GyneFix<sup>®</sup> is highly tolerated and the removal rate for pain is very low.

The removal rate for bleeding is also low with GyneFix<sup>®</sup> and is ascribed to the small surface area and the characteristics of the device already described. Clinical trials have shown that GyneFix<sup>®</sup> 200 does not increase menstrual blood loss significantly after the first few months following insertion. Although abnormal bleeding still may occur, the low removal rate for bleeding/pain is probably due to the low incidence of concomitant complaints of pain, rendering abnormal bleeding, mostly spotting, more acceptable to the woman.

Recommendation: In case of spotting or heavy bleeding occurs the use of NSAID agents may be useful.

### **3. Expulsion of GyneFix<sup>®</sup>**

With traditional IUDs, young and nulligravid/nulliparous women are particularly prone to downward migration and expulsion of the IUD. This is not the case with GyneFix<sup>®</sup>. Long-term multicenter clinical trials with GyneFix<sup>®</sup> showed lower expulsion rates both in parous and nulliparous women, as compared with expulsion rates of conventional copper IUDs, on condition that the device is properly anchored.

#### **Duration of action Safety of the anchoring system**

GyneFix<sup>®</sup> has a duration of action of 5 years.

Histology studies conducted in hysterectomy specimen up to 4 years after insertion of GyneFix<sup>®</sup> have shown that the myometrial tissue reaction at the site of the anchor is minimal (less than 1 mm), or even absent, supporting both the safety of the material and the safety of the insertion system. No transplanted endometrial tissue could be observed within the adjacent endometrium in any of these cases studied.

#### **Reversibility and Return to Fertility**

In depth studies have shown that there is no difference in conception rates between a group who had previously used a copper IUD and a group who had not used an IUD.

Fertility studies after GyneFix® removal have shown that reversal of fertility is unimpaired.

## **VII. Contraindications and caution (for IUDs in general)**

Contraindications: pregnancy, hypoplastic uterus (e.g., following use of injectable contraceptives), abnormalities of the uterus, current pelvic inflammatory disease (PID), sexual transmitted disease (STD), postpartum endometritis or infected abortion in the past 3 months, endometrial or cervical malignancy, unexplained genital bleeding, acute cervicitis, bleeding disorders (e.g., von Willebrand's disease), genital actinomycosis, known or suspected allergy to copper, diagnosed Wilson's disease, leukaemia, chronic corticosteroid or immunosuppressive therapy.

Caution: anaemia (depending on the cause), previous ectopic pregnancy, prior uterine or tubal surgery, valvular heart disease (with antibiotic cover), heavy menstrual flow, severe menstrual cramps, fibroid uterus.

## **VIII. Warnings (for IUDs in general)**

**Attention! This device is for single use only. Do not insert the device after expiry date. If re-used, the risk of pelvic infection could be increased. Do not resterilize. The device should not be used if the primary package is damaged. The device must be inserted in sterile conditions.**

### **A. Pregnancy**

#### 1. Septic abortion

Reports on traditional IUDs have indicated an increased incidence of septic abortions, associated in some instances with septicaemia, septic shock and death in patients conceiving with an IUD in place. Most of these reports have been associated with mid-trimester of pregnancy. In some cases the initial symptoms have been insidious and not easily recognized. If a pregnancy should occur with GyneFix® in situ, the device should be removed if the string is visible even though removal may increase the chances of miscarriage. If the decision is made not to remove the device or removal proves to be difficult, or if the woman elects to continue the pregnancy, she should be warned that there may be an increased risk of abortion or preterm labour and/or sepsis and she should be followed closely. It must be noted, however, that nothing of this kind has been reported with GyneFix® so far.

#### 2. Ectopic pregnancy

A woman who conceives while using an IUD should be carefully evaluated for the possibility of an ectopic pregnancy. Special attention should be directed to patients with delayed menses, heavy and prolonged bleeding and/or pelvic pain.

The risk of ectopic pregnancy with GyneFix® is lower than in women using no contraception.

### **B. Pelvic Infection**

An increased risk of PID associated with the use of IUDs has been reported. The risk appears to be greatest for nulliparous women who have a multiplicity of sexual partners. Studies conducted with GyneFix® in nulligravid and nulliparous women, however, have not shown any increased risk of PID in properly selected nulligravid and nulliparous women. These studies confirm that the non-prescription of IUDs to nulliparous women because of the fear of pelvic inflammatory disease (PID) and the potential for resulting

infertility, (a recommendation expressed throughout the 1980s), is no longer justified. These results support the current view that it is not the age or parity of the woman but her life-style which may enhance the risk of PID. If salpingitis should occur, however, the infection could result in tubal damage and occlusion, thereby threatening future fertility. Therefore, it is recommended that women be advised to look for and report symptoms of pelvic inflammatory disease. The symptoms of pelvic infection include: development of menstrual anomaly (prolonged or heavy bleeding), abnormal vaginal discharge, abdominal or pelvic pain, dyspareunia, fever. The signs and symptoms are especially significant if they occur following the first few cycles after insertion. Appropriate aerobic and anaerobic bacteriological studies should be done and antibiotic therapy initiated promptly. If the infection does not show marked clinical improvement within 24 to 48 hours, GyneFix® should be removed and treatment reassessed on the basis of the results of culture and sensitivity tests.

### **C. Perforation**

Partial or total perforation of the uterine wall or cervix has occurred with the use of conventional IUDs, most commonly during IUD insertion. Although perforations with GyneFix® are rare (approximately 0.5/1000), the possibility of perforation must be kept in mind during insertion and at the time of any subsequent examination. Visualization of the anchor following insertion should reduce the risk of perforation. If perforation occurs, the device should be removed immediately. Laparoscopy or laparotomy may be indicated to remove the device. Abdominal adhesions, intestinal penetration, obstruction, and local inflammatory reaction with abscess formation and erosion of adjacent viscera may result if the IUD is left in the abdominal cavity. These complications are extremely rare.

### **D. Effects of Copper**

Additional amounts of copper available to the body from a copper IUD may precipitate symptoms in women with undiagnosed Wilson's disease. The estimated incidence of Wilson's disease is 1 in 200,000.

### **IX. How Supplied and How Stored**

Each GyneFix® IUD is sterile (EO-sterilization) and delivered in a blister pack, together with an inserter unit and uterine sound. Store dry and between 0°C and 35°C.

### **INSERTION INSTRUCTIONS**

*(See also Video INSERTION INSTRUCTIONS on <http://www.wildemeersch.com>)*

Since GyneFix® represents a unique design in intrauterine contraception, physicians are cautioned that it is imperative for them to make themselves thoroughly familiar with the instructions for insertion before attempting placement of GyneFix®. Prior to insertion, the physician must review all aspects associated with the use of the GyneFix®. The patient should be given the opportunity to discuss fully any questions she may have. A relevant medical history should be obtained to determine conditions that might influence the selection of the IUD as a method of contraception. Physical examination should include a pelvic examination and, if indicated, a 'Pap' smear, and appropriate tests for other forms of genital disease. Pregnancy should be ruled out prior to insertion. If the woman is anxious, the use of local intra-cervical or local/ regional anaesthesia should be considered. The use of misoprostol prior to fitting of the IUD may sometimes also be

useful to dilate the cervical canal. IUD providers should be aware that the popularity of the IUD could be much improved if attention is given to pain relief during IUD insertion.

**GyneFix® is a frameless, flexible IUD, which is fixed in the fundus of the uterus.**

**The following are some pertinent recommendations:**

- The thickness of the fundus should be at least 1.1 cm. In case of doubt the provider should evaluate the fundal thickness by ultrasound.
- The use of a forceps ,18 or 19 cm in length, Allis or Pozzi forceps, is recommended as this facilitates insertion.
- Trim the thread with a scalpel.
- Attention to comfort during insertion is very important. If necessary, perform cervical priming by, for example, Cytotec®: 2 tablets of 200mg orally, 3 hours before the insertion. In addition, you can also use intra-cervical anesthesia, preferably with a dental syringe.
- Following delivery, especially in breastfeeding mothers, the insertion should be postponed until three months after birth.
- Women using the 3-monthly injectable should not use GyneFix® because of atrophy of the uterus. These women can only use GyneFix® when the uterus has regained a normal size.
- After the insertion, the woman should not have intercourse, not use tampons nor Mooncup within 14 days.

If properly inserted, GyneFix® has several important advantages when compared to traditional IUDs. It is therefore essential that doctors become familiar with this new insertion technique. This leaflet explains how proficiency can be acquired in a minimum of time.

**Remark:**

- The physician is responsible to inform the patient about possible side-effects.
- Never insert GyneFix® without proper training.
- After you have participated in a training course, use HUT® (the “Home Uterine Trainer”) to practice at home before you start inserting GyneFix® in your patients. The instructions below concern the insertion procedure in HUT®. The procedure is identical with the actual insertion in your patients.

**NOTE:**

Experience has shown that most mistakes are made by not following correctly the instructions (**see main instructions in bold**).

It is very important to measure the fundal thickness by ultrasound. The thickness should be at least 1.1 cm.

**Step 1:** Start with a gynaecological examination and disinfection of the vagina. Place the Allis or Pozzi forceps on the cervix (horizontally or vertically as you prefer).

**Step 2:** Put your hand on the forceps (between the first and second finger) as shown on the figure. **Make sure your thumb is free.**

**Step 3:** **Sound the uterus** with the sound delivered in the package: it is very important to know depth and direction of the uterine cavity. Successful sounding of the uterus also provides confidence: if you can make contact with the fundus with the sound, you will also be able to make contact with the GyneFix® applicator.

**Step 4: Grasp the applicator as shown on the figure.** Take the applicator with the other hand than the one you are holding the Allis Forceps. Take the applicator between your thumb and index finger out of the package. You can hold the applicator correctly, make sure the needle is at the same level of the insertion. **Maintain this grip until you have reached the fundus (See step 7).**

**Step 5:** Insert the applicator in the uterus **until it touches the fundal wall** of the uterus. The depth should be identical as the sound length. Do not change the adjustment ring on the application tube. It only serves to fix the applicator in the package, and remains at 10 cm.

**Step 6:** Immediately following contact with the fundus, **put your thumb of the hand holding the forceps at the end of the handle** (see figure) to keep contact with the fundus before anchoring. Don't push just keep contact.

**Step 7:** Remove your hand from the applicator now.

**Step 8: Position your thumb and index finger on the handle as shown.**

**Step 9:** Focus now on the distance between the handle and the end of the insertion tube (see figure).

**Step 10:** Move the handle forward, feeling the anchoring knot penetrating the fundus. Full penetration is obtained when the handle touches the tube. **It is not necessary to push hard or to push a second time.**

**Step 11:** When practising on HUT<sup>®</sup>, remove the tail which is fixed in the slot on the handle as shown (note: cut the tail with a scalpel when you insert GyneFix<sup>®</sup> in your patients). **Make sure that your thumb is still positioned at the end of the handle to avoid traction on the anchor.**

**Step 12:** Remove your thumb now and withdraw the handle from the insertion tube.

**Step 13:** Withdraw the tube slowly while rotating it.

**Step 14:** Optional: Stroke the tail softly to feel the attachment. Do not pull the tail immediately after the procedure. Trim the tail at 3-4cm.

**Step 15:** The presence of the tiny metal tube allows visualization on ultrasound. This examination provides important information on the proper placement of the anchor. After insertion, look for the visualization marker and measure the distance to the serosa to check if the insertion was done correctly. In case the insertion was not done correctly, remove the device immediately

## **REMOVAL**

GyneFix<sup>®</sup> is removed simply by a short pull at the tail (Note: asking the patient to cough reduces pain sensations).

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